## **Pipette Shakers**

LW Scientific pipette shakers easily achieve consistent results in specimen separation and are especially effective for blood preparation. The LWS 6 place shaker has the option of either continuous motion or timed motion using the 15 minute timer. Compact size and durability make the LW Scientific pipette shaker convenient for any laboratory.

## SPECIFICATIONS

Portable and compact Capacity – 6 places Holds 5ml - 7ml pipettes Fixed speed - 2500 rpm 15-minute timer Suction-cupped feet Space-saving design LED indicator light Durable construction 110v or 220v One year warranty



Model #

SHL-PPF7-06F1

Pipette shaker 6 places



P 770.270.1394 F 770.270.2389

865 Marathon Parkway Lawrenceville GA 30046

LWScientific.com